


Insta-TEST[®]

25 Test Strips

IRON

YaliTech[®]
INSTRUMENTS
FOOD & WATER

LaMotte
Just DIP & READ!



Immerse strip for 2 seconds.
Remove with pad face up.
SHAKE ONCE TO REMOVE EXCESS WATER. Wait 60 seconds and compare to color chart.



View from above.

Iron ppm



LaMotte Company,
Chestertown, MD 21620 USA
www.lamotte.com

Code 2935-G



- Keep wet fingers out of the bottle.
- Close cap tightly after removing strip.
- In pools, immerse strip to a depth of 6".
- In spas, turn off jets and immerse strip to a depth of 6".
- For best results read in natural daylight.
- Store in a cool, dry place.



6 37395 46617 5

Insta-TEST®

25 Test Strips
25 Languettes d'analyse
25 Tiras de análisis

- COPPER
- CUIVRE
- COBRE

 LaMotte

Immerse strip for **2 seconds**. Remove strip and compare to color chart. **REMOVE EXCESS WATER. Wait 30 seconds** and compare to color chart.

- ! Keep wet fingers out of the bottle.
- ! Close cap tightly after removing strip.
- ! In pools or spas, immerse strip to a depth of 6" (15 cm).
- ! For best results read in natural daylight.
- ! Store in a cool, dry place.

Plongez la bande dans l'eau pendant **2 secondes**. Sortez-la de manière à retirer l'excès d'eau. **ATTENDEZ 30 SECONDES POUR ÉLIMINER L'EXCÈS D'EAU**. Attendez **30 secondes** et comparez avec la carte de couleurs.

- ! Maintenez les doigts humides hors de la bouteille.
- ! Bien refermer le bouchon après avoir retiré la bande réactive.
- ! Pour les piscines et les spas, plonger la bande à une profondeur de 15 cm (6").
- ! Pour un meilleur résultat, réaliser la lecture à la lumière naturelle.
- ! Conserver dans un endroit frais et sec.

Copper, Cuivre, Cobre ppm



0 0.3 0.6 1.0 3.0

- ! **Sumerja la tira durante 2 segundos**. Sáquela de modo que la almohadilla quede hacia arriba. **SACUDALA UNA VEZ PARA ELIMINAR EL EXCESO DE AGUA**. Espere **30 segundos** y compare con la carta de color.
- ! Mantenga los dedos húmedos fuera de la botella.
- ! Cierre bien el tapón después de sacar la tira reactiva.
- ! En piscinas o spas, sumerja la tira a una profundidad de 15 cm (6").
- ! Para un mejor resultado, realice la lectura con luz natural.
- ! Conserve en un lugar fresco y seco.

Code 2991 - G
802 Washington Ave.
Chesapeake, MD 21620 USA
410-778-3100
www.lamotte.com/insta



Drop holding box and dotted tolerance line